

Assay background

This protocol is designed to evaluate the efficiency of peptide-HLA-I interaction and complex formation. The assay is based on detecting the β_2 -microglobulin (β_2m) light chain subunit of recombinant HLA class I (HLA-I) complexes, where the heavy chain has been biotin tagged. These tagged complexes are subsequently captured by streptavidin coated beads, labelled with PE-conjugated anti-human β_2m , and analyzed by flow cytometry. Since peptide-HLA-I complex formation is entirely peptide dependent, bead-associated signals will only be detected if the peptide in question supports the folding of the HLA-I allotype of interest; peptides that efficiently support folding will give strong signals whereas peptides that support folding sub-optimally, or not at all, will give moderate to non-detectable signals.

Materials and equipment, you need to provide

Dilution buffer: PBS with 5% glycerol.

FACS buffer: PBS with 1%BSA (or FCS) and 0.01% NaN₃.

Streptavidin coated beads: (6-8 μ m beads from Spherotech Cat.: SVP-60-5)

BBM.1-PE: (anti-human β_2m) (Santa Cruz Cat# sc-13565 PE)

Sealing Tape: (Thermo Scientific Nunc™ cat.: 236366)

96-well plate: U-bottom shape

Centrifuge with a plate rotor


Rocking table

Flow cytometer

Protocol

- The folded HLA complexes are prepared according to protocol: "HLA-I tetramer production". The folding setup should include a **positive control** (provided with the easYmer®): a peptide that is known to support folding of the HLA molecule of interest, and a **negative control**: without peptide.
- In the easYmer® folding set-up, the highest achievable concentration of the folded complex is 500nM.
- After completed folding incubation the complex formation can be evaluated.
- Prepare sufficient dilution buffer (PBS, 5% glycerol) for the whole assay.

- Dilute each of the folded complexes to give 75 μ L of a 40nM solution (e.g. for a 500nM complex: 6 μ L folded complex in 69 μ L dilution buffer).
- For all samples and **positive and negative controls**, transfer 50 μ L of this pre-dilution (prepared in 5.) to the first tube. Make three subsequent serial 3-fold dilutions (50 μ L in 100 μ L dilution buffer), according to the figure below.



- Transfer 40 μ L of each these dilutions to the wells in a U-bottom shape 96-well plate, as suggested below. Also prepare a background well (BLANK): transfer 40 μ L of dilution buffer.
- Prepare a sufficient volume of a 45-fold dilution of the streptavidin coated beads (6-8 μ m; SVP-60-5) in dilution buffer. Transfer 20 μ L of the diluted bead suspension to each well.

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLANK		P-1		S1-1		S3-1		S5-1		S7-1	
B			P-2		S1-2		S3-2		S5-2		S7-2	
C			P-3		S1-3		S3-3		S5-3		S7-3	
D												
E			N-1		S2-1		S4-1		S6-1		S8-1	
F			N-2		S2-2		S4-2		S6-2		S8-2	
G			N-3		S2-3		S4-3		S6-3		S8-3	
H												

BLANK : No complex


P1-3 : Positive control dilutions (HLA with know peptide)

N1-3 : Negative control dilutions (HLA without peptide)

S1-S8 : Sample dilutions (complexes to evaluate)

- Mix well and seal the plates with Sealing Tape to avoid well to well contamination.
- Incubate the plate on a rocking table at 37°C for 1h.
- Remove the Sealing Tape and wash by adding 160µl FACS buffer.
- Spin the plate at 700g for 3min and flip out the supernatant.
- Resuspend the beads in 200µl FACS buffer.
- Spin the plate at 700g for 3min and flip out the supernatant.
- Wash two more times by repeating step 13. and 14.
- During the above washing steps, prepare a 200-fold dilution of the PE-labeled anti-human β_2m monoclonal antibody BBM.1 in FACS buffer
- Resuspend the beads in 50µL antibody solution per well.
- Incubate the plate for 30 min. at 4°C.
- Wash by adding 150µl FACS buffer. Spin the plate at 700g for 3min and flip out the supernatant.
- Resuspend the beads in 200µl FACS buffer. Spin the plate at 700g for 3min and flip out the supernatant.
- Wash two more times by repeating step 19. and 20.
- Resuspend the beads in 200µl FACS buffer, and analyze on a Flow cytometer.

Example of the Flow cytometry-based assay:



Flow cytometry-based detection of 4 different peptide-HLA-A*02:01 complexes.

Complexes of A*02:01 and 4 different peptides, and a negative control (No Peptide), were folded. CMV pp65₄₉₅₋₅₀₃ (NLVPMVATV) a known HLA-A*02:01 restricted epitope was used as positive control. The three other peptides are based on their A*02:01 binding stability categorized as good binder ($T_{1/2}$ 6.5h), intermediate binder ($T_{1/2}$ 3.5h), and low binder ($T_{1/2}$ 0.7h). Three dilutions of the folded complexes were analysed in the flow cytometry-based assay. The X-axis gives the complex concentration if complete folding is achieved.