

Assay background

This protocol is designed to evaluate the efficiency of peptide-HLA-I interaction and complex formation. The assay is based on detecting the β_2 -microglobulin (β_2m) light chain subunit of recombinant HLA class I (HLA-I) complexes, where the heavy chain has been biotin tagged. These tagged complexes are subsequently captured by streptavidin coated beads, labelled with PE-conjugated anti-human β_2m , and analyzed by flow cytometry. Since peptide-HLA-I complex formation is entirely peptide dependent, bead-associated signals will only be detected if the peptide in question supports complex formation of the HLA-I allotype of interest; peptides that efficiently support complex formation will give strong signals whereas peptides that support sub-optimal complex formation, or not at all, will give moderate to non-detectable signals.

Materials and equipment, you need to provide

Dilution buffer: PBS with 5% glycerol.

FACS buffer: PBS with 1%BSA (or FCS) and 0.01% NaN_3 .

Streptavidin coated beads: (6-8 μ m beads from Spherotech Cat.: SVP-60-5)

BBM.1-PE: (anti-human β_2m) (Santa Cruz Cat# sc-13565 PE)

Sealing Tape: (Thermo Scientific Nunc™ cat.: 236366)

96-well plate: U-bottom shape

Centrifuge with a plate rotor


Rocking table

Flow cytometer

Protocol

- The peptide loaded HLA complexes are prepared according to protocol: "HLA-I tetramer production". The setup should include a **positive control** (provided with the easYmer®): a peptide that is known to support complex formation of the HLA molecule of interest, and a **negative control**: without peptide.
- In the easYmer® set-up protocol, the highest achievable concentration of the complex is 500nM.
- Prepare sufficient dilution buffer (PBS, 5% glycerol) for the whole assay.
- Dilute each complex to give 75 μ L of a 40nM solution (e.g. for a 500nM complex: 6 μ L complex in 69 μ L dilution buffer.

- For all samples and **positive and negative controls**, transfer 50 μ L of this pre-dilution (prepared in 5.) to the first tube. Make three subsequent serial 3-fold dilutions (50 μ L in 100 μ L dilution buffer), according to the figure below.



- Transfer 40 μ L of each these dilutions to the wells in a U-bottom shape 96-well plate, as suggested below. Also prepare a background well (BLANK): transfer 40 μ L of dilution buffer.
- Prepare a sufficient volume of a 45-fold dilution of the streptavidin coated beads (6-8 μ m; SVP-60-5) in dilution buffer. Transfer 20 μ L of the diluted bead suspension to each well.
- Mix well and seal the plates with Sealing Tape to avoid well to well contamination.
- Incubate the plate on a rocking table at 37°C for 1h.

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLANK		P-1		S1-1		S3-1		S5-1		S7-1	
B			P-2		S1-2		S3-2		S5-2		S7-2	
C			P-3		S1-3		S3-3		S5-3		S7-3	
D												
E			N-1		S2-1		S4-1		S6-1		S8-1	
F			N-2		S2-2		S4-2		S6-2		S8-2	
G			N-3		S2-3		S4-3		S6-3		S8-3	
H												

BLANK : No complex


P1-3 : Positive control dilutions (HLA with know peptide)

N1-3 : Negative control dilutions (HLA without peptide)

S1-S8 : Sample dilutions (complexes to evaluate)

10. Remove the Sealing Tape and wash by adding 160µl FACS buffer.
11. Spin the plate at 700g for 3min and flip out the supernatant.
12. Resuspend the beads in 200µl FACS buffer.
13. Spin the plate at 700g for 3min and flip out the supernatant.
14. Wash two more times by repeating step 13. and 14.
15. During the above washing steps, prepare a 200-fold dilution of the PE-labeled anti-human β_2m monoclonal antibody BBM.1 in FACS buffer
16. Resuspend the beads in 50µL antibody solution per well.
17. Incubate the plate for 30 min. at 4°C.
18. Wash by adding 150µl FACS buffer. Spin the plate at 700g for 3min and flip out the supernatant.
19. Resuspend the beads in 200µl FACS buffer. Spin the plate at 700g for 3min and flip out the supernatant.
20. Wash two more times by repeating step 19. and 20.
21. Resuspend the beads in 200µl FACS buffer, and analyze on a Flow cytometer.

Example of the Flow cytometry-based assay:



Flow cytometry-based detection of 4 different peptide-HLA-A*02:01 complexes.

Complexes of A*02:01 and 4 different peptides, and a negative control (No Peptide), were prepared. CMV pp65₄₉₅₋₅₀₃ (NLVPMVATV), a known HLA-A*02:01 restricted epitope, was used as positive control. The three other peptides are based on their A*02:01 binding stability categorized as good binder ($T_{1/2}$ 6.5h), intermediate binder ($T_{1/2}$ 3.5h), and low binder ($T_{1/2}$ 0.7h). Three dilutions of the complexes were analysed in the flow cytometry-based assay. The X-axis gives the complex concentration if complete loading is achieved.