

HLA-C*03:04 easYmer kit

| | |
|-----------------------|---|
| Catalog# | 1117-01 |
| Size | 20, 50, 150, or 500 tests |
| easYmer | HLA-C*03:04 - peptide receptive, biotinylated in Tris/Maleate pH 7 with 30% Glycerol |
| Folding buffer | Tris/Maleate pH 7 |
| Peptide | RAVEPGTVL, a stable HLA-C*03:04 binder Positive control for evaluation analysis of peptide-HLA folding |
| Application | easYmers can be used to generate monomers with your choice of peptide. The monomers can easily be tetramerized with fluorophore conjugated streptavidin and used to analyse T cells by flow cytometry. The easYmer reagent can also be used to evaluate specific peptide-HLA I binding. |

For Research Use Only (RUO)